

3-8% Tris-Acetate Midi Protein Gels, 26-well, w/adapters

Product Information

Product Name

3-8% Tris-Acetate Midi Protein Gels, 26-well, w/adapters

Category

SDS-PAGE

Research

Tris-Acetate

Function

Denaturing, Native

Unit Size

10 gels (1 box)

Sample Loading Volume

up to 15 μ L

Mode of Separation

Molecular Weight

Shipping Condition

Wet Ice

Gel Thickness

1.0 mm

Well Design

1D Well

Length (Metric)

8 cm

Length (Metric) Gel

3-8% Tris-Acetate Midi Protein Gels, 26-well, w/adapters

8 cm

Shelf Life

6 Months

Width (Metric)

13 cm

Width (Metric) Gel

8 cm

Product Type

Gradient Page

Gel Size

Midi

Separation Range

36 to 500 kd

Gel Percentage

3 to 8%

Gel Type

Tris-Acetate

For Use With (Equipment)

Midi Gel Tank

Separation of

Protein

Wells

26-well, w/adapters

Quantity

3-8% Tris-Acetate Midi Protein Gels, 26-well, w/adapters

1 Box

Contents & Storage

One box contains 10 gels. Store in refrigerator (2–8°C). Do not freeze. Shelf life is 6 months.
