

4-12% Tris-Glycine Plus Midi Gels 26-well, w/adapters

Product Information

Product Name

4-12% Tris-Glycine Plus Midi Gels 26-well, w/adapters

Category

SDS-PAGE

Research

Tris-Glycine

Function

Denaturing, Native

Unit Size

10 gels (1 box)

Shipping Condition

Wet Ice

Gel Thickness

1.0 mm

Shelf Life

12 Months

Gel Size

Midi

Separation Range

20 to 250 kd

Gel Percentage

4 to 12%

Gel Type

4-12% Tris-Glycine Plus Midi Gels 26-well, w/adapters

Tris-Glycine

For Use With (Equipment)

XCell4 SureLock Midi

Separation of

Protein

Wells

26-well, w/adapters

Quantity

10/Pk.

Contents & Storage

One box containing 10 gels and 10 adapters. Store in refrigerator (2–8°C). Do not freeze.
