

## 8-16% Tris-Glycine Plus Midi Gels 26-well, w/adapters

### Product Information

---

**Product Name**

8-16% Tris-Glycine Plus Midi Gels 26-well, w/adapters

---

**Category**

SDS-PAGE

---

**Research**

Tris-Glycine

---

**Function**

Denaturing, Native

---

**Unit Size**

10 gels (1 box)

---

**Shipping Condition**

Wet Ice

---

**Gel Thickness**

1.0 mm

---

**Shelf Life**

12 Months

---

**Gel Size**

Midi

---

**Separation Range**

6 to 250 kd

---

**Gel Percentage**

8 to 16%

---

**Gel Type**

---

## 8-16% Tris-Glycine Plus Midi Gels 26-well, w/adapters

Tris-Glycine

---

**For Use With (Equipment)**

XCell4 SureLock Midi

---

**Separation of**

Protein

---

**Wells**

26-well, w/adapters

---

**Quantity**

10/Pk.

---

**Contents & Storage**

One box containing 10 gels and 10 adapters. Store in refrigerator (2–8°C). Do not freeze.

---