

8-16% Tris-Glycine Plus Midi Gels, 12+2-well, w/adapters

Product Information

Product Name

8-16% Tris-Glycine Plus Midi Gels, 12+2-well, w/adapters

Category

SDS-PAGE

Research

Tris-Glycine

Function

Denaturing, Native

Unit Size

10 gels (1 box)

Shipping Condition

Wet Ice

Gel Thickness

1.0 mm

Shelf Life

12 Months

Gel Size

Midi

Separation Range

6 to 250 kd

Gel Percentage

8 to 16%

Gel Type

8-16% Tris-Glycine Plus Midi Gels, 12+2-well, w/adapters

Tris-Glycine

For Use With (Equipment)

XCell4 SureLock Midi

Separation of

Protein

Wells

12+2-well, w/adapters

Quantity

10/Pk.

Contents & Storage

One box containing 10 gels and 10 adapters. Store in refrigerator (2–8°C). Do not freeze.
