

Anti-HA Agarose

Product Information

Product Name	
Anti-HA Agarose	
Category	
Protein Purification Kit	
Purification Target	
HA-tagged Proteins	
Unit Size	
5 mL	
Column Type	
Agarose Resin, Affinity	
Stationary Phase	
Antibody	
Format	
Bottle(s)	
Contents & Storage	
5 mL of settle resin, provided as 10 mL of 50% slurry. Store at 2–8°C.	

5 mL of settle resin, provided as 10 mL of 50% slurry. Store at 2–89

Tel: 1-631-559-9269 1-516-512-3133

Fax:1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA